

Alkanes-Mix - 34 components: 1000ug/ml each in n-Hexane (5x1ml ampules)

Chromatogram and analytical data included in the COA

CRM, ISO 17034 and ISO 17025

d

34 components: 1000mg/l each of

n-Heptane (C7) [CAS:142-82-5]
n-Octane (C8) [CAS:111-65-9]
n-Nonane (C9) [CAS:111-84-2]
n-Decane (C10) [CAS:124-18-5]
n-Undecane (C11) [CAS:1120-21-4]
n-Dodecane (C12) [CAS:112-40-3]
n-Tridecane (C13) [CAS:629-50-5]
n-Tetradecane (C14) [CAS:629-59-4]
n-Pentadecane (C15) [CAS:629-62-9]
n-Hexadecane (C16) [CAS:544-76-3]
n-Heptadecane (C17) [CAS:629-78-7]
n-Octadecane (C18) [CAS:593-45-3]
n-Nonadecane (C19) [CAS:629-92-5]
n-Eicosane (C20) [CAS:112-95-8] ;
n-Heneicosane (C21) [CAS:629-94-7]
n-Docosane (C22) [CAS:629-97-0]
n-Tricosane (C23) [CAS:638-67-5]
n-Tetracosane (C24) [CAS:646-31-1]
n-Pentacosane (C25) [CAS:629-99-2]
n-Hexacosane (C26) [CAS:630-01-3]
n-Heptacosane (C27) [CAS:593-49-7]
n-Octacosane (C28) [CAS:630-02-4]
n-Nonacosane (C29) [CAS:630-03-5]
n-Triacontane (C30) [CAS:638-68-6]
n-Hentriacontane (C31) [CAS:630-04-6]
n-Dotriacontane (C32) [CAS:544-85-4]
n-Tritriacontane (C33) [CAS:630-05-7]
n-Tetratriacontane (C34) [CAS:14167-59-0]
n-Pentatriacontane (C35) [CAS:630-07-9]
n-Hexatriacontane (C36) [CAS:630-06-8]
n-Heptatriacontane (C37) [CAS:7194-84-5]
n-Octatriacontane (C38) [CAS:7194-85-6]
n-Nonatriacontane (C39) [CAS:7194-86-7]
n-Tetracontane (C40) [CAS:4181-95-7]
in n-Hexane

Conservation 12 months.

Store in a refrigerator at temperatures between 2°C to 8°C